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
1. Load an existing isotherm | 2. Display of expected isotherm | 3. Input of measurement range | 4. Automatic optimization of conditions such as gas introduction volume |

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

1. Room temperature change | 2.  $P_{ref}$  | 3.  $P_{smp}$  | 4. Initial pressure | 5. The sample tube temperature changes due to the change of  $LN_2$  level | 6. Liquid nitrogen temperature changes due to oxygen dissolution | 7. AFSM | 8. Sample

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1. Reference cell for measuring fluctuations within free space continuously | 2. Reference cell | 3. Sample cell | 4. Liquid level 1 | 5. Liquid level 2

AFSM (He-based method)

1. **Sample exhaust**
2. Free space (@room temperature) measurement with sample
3. Sample exhaust
4. **Reference tube,  $p_0$  tube gas introduction**
5. Free space (@adsorption temperature) measurement with sample
6. Sample exhaust
7. **Adsorption measurement, free space change acquisition**
8. **Sample exhaust**
9. **Purge**

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The versatility of BELSORP equipment is truly world leading. The numerous features and capabilities are complemented by BELCONTROL the intuitive and user-friendly operation software. It guides the user step-by-step through the analysis process. This includes the setup of analysis conditions, executing the measurements, when to fill and setup the liquid nitrogen or other bath, when to replace the gas cylinder, the degassing steps, and much more. The software is designed to make the instrument accessible and operable to everyone, including inexperienced users.

For inexperienced users or for measurements of unknown samples, BELCONTROL only requires basic sample information (name, mass, etc.), pre-treatment conditions (if not performed externally) and the measurement range.

Detailed control of the configuration and measurement settings is possible to optimize the measurement conditions (e.g. dosing settings, equilibrium criteria, leak test option, etc.). This allows the user to fully customize the sample analysis to his needs.

---

**BELCONTROL**

Quick BET	igen	Multi-point BET surface area in less than 20 minutes
Helium-Free Measurement	igen	AFSM <sup>TM</sup> 2 enables He-free measurements with unmatched accuracy
Adsorption Kinetics	opcionális	Rate of adsorption measurements for diffusion analysis

## MORE BELCONTROL FEATURES

- | Overlaying adsorption / desorption isotherms and comparing the measured data between the various ports during measurement
- | All pressures, temperatures, valve actuations, etc. are stored in trend data, allowing for immediate examination
- | A system check function is available for diagnosing the instrument status
- | E-mail notification automatically transmits measurement status and results
- | Interactive program in Japanese or English ensures easy, reliable operation
- | Extensive help functions, including step-by-step instructions during operation

BET SPECIFIC SURFACE AREA & PORE SIZE ANALYZER BELSORP MINI X

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Microtrac offers various accessories for its entire gas & vapor adsorption product line.

### STANDARD CONSUMABLE GOODS

The standard consumables comprise sample cells, filler rods, filters, O-rings, caps and weighing platforms. Various sizes of sample cells, quick seals, and other optional consumables are also available.

### WATER BATH

Easily swap the Dewar with a water bath for measurement temperatures ranging from  $-10^{\circ}\text{C}$  to  $70^{\circ}\text{C}$ . A refrigerated / heated circulator is required.

## HEATER

In-situ pretreatment of the sample from 50°C to 450°C. With this option, samples can be pretreated directly at the measurement port, without the need to transfer the sample cell before analysis. Particularly useful for sensitive samples.

## GAS SELECTOR


The gas selector enables the connection of up to 4 adsorptives simultaneously (1 x Helium + 4 x adsorptive).

BET SPECIFIC SURFACE AREA & PORE SIZE ANALYZER BELSORP MINI X

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Microtrac's gas adsorption analyzers are used in a variety of fields. These include catalysts, batteries, fibers, polymer materials, zeolite, fuel cells, chemicals, pigments, cosmetics, MOF / PCP, magnetic powders, separation membranes, filters, toners, cement, ceramics, semi-conductors, and many more.

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- | Szilika
- | MOF / PCP
- | pigmentek
- | kozmetikumok

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The principles of volumetric adsorption measurement are explained in our knowledge base:

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


When developing the software, we focused on simplifying the operations and equipped it with many functions to increase labor productivity. Since the BELSORP instruments offer so many features, it is important to make those features easy to use. Our software will guide you step-by-step for the implementation of several procedures, i.e. execution of measurements, replacement of gas cylinders and purging of the manifold. The step-by-step instructions allows for users of all levels of experience to use the machinery safely and effectively. Depending on the level of user-experience two modes are offered, Simple Mode and Professional Mode. Simple Mode enables an easy operation by decreasing the customizability of the measurement – the only information you will need is the sample labeling, selecting pretreatment conditions (skippable if externally done) and measurement points. This low maintenance measurement is ideal for measuring unknown samples. If a prior measurement with comparable sorption behavior is available, the GDO function can be used to reduce the measurement time. Professional Mode offers the user a detailed configuration of the measurement, i.e. dosing settings, equilibrium criteria, leak-check option, and more.

## GAS ADSORPTION MEASUREMENTS OF POROUS & NON-POROUS MATERIALS: MOFS, ZEOLITES, CARBONS, AND MORE

### EVALUATION SOFTWARE BELMASTER

A sorption isotherm is defined as the relationship between the adsorbed amount onto an adsorbent and the equilibrium pressure of the adsorptive gas – commonly related to the saturation vapor pressure – at constant temperature. The gas sorption isotherm (e.g. nitrogen) delivers information about the specific surface area, pore size distribution and pore volume of measured material. In the following graph some exemplary sorption isotherms are shown.

The specific surface area (SSA) refers to the accessible surface area of the sample and is of great importance in adsorption, homogeneous and heterogeneous catalytic reactions. The specific surface area can be calculated according to ISO 9277 with the BET method (BET: Brunauer, Emmett and Teller) or Langmuir method. The following graph shows an exemplary evaluation of the specific surface area according to the BET method in our BELMASTER software:

Select the right pressure range (multi point BET) or a measurement point (single point BET) and the surface area will be calculated automatically. Further, our BELMASTER software also provides the BET surface area calculation according to ISO 9277 Annex C (also known as Rouquerol-plot), which is recommended for microporous materials.

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## C KONSTANS

### AFSM™


The following measurements are conducted using our patented, Advanced Free Space Measurement (AFSM™). This method allows for high reproducibility by using a reference cell to follow the changes in the free space volume. Our instrument measures the initial free space of the sample cell and reference cell. Since the change of the free space is equal in both cells (same environmental conditions), the change can be continuously followed by the pressure change of the reference cell.

In a product test, a certified reference material CRM-170 (nonporous alpha alumina; certified value of nitrogen BET specific surface area SSA:  $1.05 \pm 0.05 \text{ m}^2/\text{g}$ ) is used to check the accuracy of measurement in low total surface area (TSA) measurements. Although the TSA is reduced from approx.  $2 \text{ m}^2$  to lower than  $0.4 \text{ m}^2$ , the determined BET specific surface areas SSA stay almost identical.

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
Theory	Interaction Surface and Gas	Adsorbate	Applicable Pore Size Range
BJH, CI, DH, INNES method	Kelvin equation (Surface tension and contact angle)	Bulk liquid density	> 2 nm Meso and macropore
HK, SF, CY method	Lennard-Jones potential (Interaction and repulsion force)	Bulk liquid density	0.4 - 2 nm Micropore
NLDFT, GCMC	Statistical thermodynamic model		0.35 - 500 nm Whole pore range

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
In recent years, attention has been focused on pore structure evaluation methods using computer simulations, such as the novel pore distribution analysis NLDFT (Non-localized Density Functional Theory) and GCMC (Grand Canonical Monte Carlo) method, which can analyze micropores to meso- and macropores using a unified theory. Pore size distributions obtained from the same adsorption isotherm are different between classical and novel PSD analyses, and even in between novel methods because the filling pressure obtained from each theory is different. Microtrac provides evaluation methods which cover a wide range of pore sizes and various adsorbates, such as N<sub>2</sub> (77.4 K), Ar (87.3 K), and CO<sub>2</sub> (298 K). It uses NLDFT / GCMC kernels of slit, cylinder, and cage pore models with carbon and metal oxide surface atoms, resulting in the most appropriate description of porous materials.

The BELMASTER software allows easy comparison between experimental and simulated isotherms, with the simulated isotherm serving as a basis for the PSD calculation.

In the following part an example for pore size distribution calculation by BJH method is given:

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